	<h1>WEBTank</h1> <p>Manuale d'USO</p> <p>Codifica: MDU_WT-0519</p>	Ed.01	Rev.04	06/03/21
---	--	-------	--------	----------

WEBTank

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1 - [Product = THERMOVAC TTR 91 N - 2SP]; [Part No. = 230040V02]

Technical Data		THERMOVAC Transmitter		
		TTR 91 N	TTR 91 N S	TTR 96 N S
Measurement range	mbar (Torr)	5.0 x 10 ⁻⁴ to Atmosphere (3.75 x 10 ⁻⁴ to Atmosphere)		
Measurement uncertainty of reading (typical) % ¹	mbar	5 x 10 ⁻⁴ to 1 x 10 ⁻³ ±10 % 1 x 10 ⁻³ to 100 ±5 % 100 to atm ±25 %		
Repeatability of reading (typical) % ¹	mbar	1 x 10 ⁻³ to 100 ±2 %		
Sensor		MEMS-Pirani	MEMS-Pirani	Coated MEMS-Pirani
Measurement principle		Thermal conductivity according to Pirani		
Supply voltage	V DC	9 to 30		
Power consumption	W	< 1.2		
Electrical connection	V	FCC 68, RJ45 (analog) / Sub-D 15 PIN (digital)		
Analog output	V DC	V _{out} = log10(P _{meas}) x 1.288 + 6.143		
Resolution	bit	0.61 to 10		
Impedance	Ω	100		
Update rate	Hz	16		
Interfaces		FCC 68, RJ45 (analog) / RS232 (digital)		
Set point		LED-ring (160°)		
Range	mbar (Torr)	2,7 x 10 ⁻⁴ to 1000 (2,0 x 10 ⁻⁴ to 750)	2,7 x 10 ⁻⁴ to 1000 (2,0 x 10 ⁻⁴ to 750)	2,7 x 10 ⁻⁴ to 1000 (2,0 x 10 ⁻⁴ to 750)
Relay		0	2	2
Relay contact rating		1 A at 30 V AC / DC, resistive load	1 A at 30 V AC / DC, resistive load	1 A at 30 V AC / DC, resistive load
Relay contact resistance, max.	mΩ	100	100	100
Relay contact endurance, min.		100 000	100 000	100 000
1.0 A at 30 V DC load		2 000 000	2 000 000	2 000 000
0.2 A at 30 V DC load				
Status indicators		LED-ring (160°)		
Max. cable length	m	100		
Overpressure limit (abs.)	bar	6 (10 bar for inert gases)		
Operating temperature range ²	°C (°F)	0 to 60 (32 to 140)		
Storage temperature range	°C (°F)	-20 to +85 (-4 to 149)		
Max. bakeout temperature	°C (°F)	85 (185), non-operating (bakeable version up to 250 (482))		
Max. rel. humidity	% r.c.	≤ 95		
Installation orientation		Any		
Materials exposed to vacuum		304 stainless steel, Tin, Gold, Vitor®	304 stainless steel, Tin, Gold, Vitor®	304 stainless steel, Vitor®, Polyene-HT®
Dead volume (DN 16 ISO-KF)	cm ³	2.80		
Weight (DN 16 ISO-KF)	g	170		
Protection class	IP	40		
CE certification		EMC Directive 2014/53/EEC		
Controller type		DISPLAY ONE / TWO / THREE and GRAPHIC ONE / TWO / THREE		

¹ Accuracy and repeatability are typical values measured with Nitrogen gas at ambient temperature after zero adjustment
² There may be minimal deviation tolerances in the range of 40 - 85 °C

2 - [Product = THERMOVAC TTR 101]; [Part No. = 230350V01]

WEBTank

Manuale d'USO

Codifica: MDU_WT-0519

Ed.01

Rev.04

06/03/21

Technical Data

THERMOVAC Transmitter TTR 101 / TTR 101 S2 / TTR 101 S2 PB

Measurement principle 10 to 1500 mbar (7.5 to 1125 Torr) 1 to 10 mbar (0.75 to 7.5 Torr) 5 x 10 ⁻⁴ to 1 mbar (3.75 x 10 ⁻⁴ to 0.75 Torr)	Capacitive diaphragm sensor Mixed range Thermal conductivity acc. to Pirani
Display range mbar (Torr)	5 x 10 ⁻⁵ to 1500 (3.8 x 10 ⁻⁵ to 1125)
Measurement range mbar (Torr)	5 x 10 ⁻⁴ to 1500 (3.8 x 10 ⁻⁴ to 1125)
Measurement uncertainty (N ₂) 5 x 10 ⁻⁴ to 1 x 10 ⁻³ mbar 1 x 10 ⁻³ to 100 mbar 100 to 950 mbar 950 to 1050 mbar	± 50% of reading ± 15% of reading ± 5% of reading ± 2,5% of reading
Repeatability (N ₂) 1 x 10 ⁻³ to 1100 mbar	± 2% of reading
Output signal (measurement signal) Voltage range Measurement range Error signal Relationship between voltage and pressure	0 to 10.23 V +0.61 to +10.23 V 0 V (factory default) 1.286 V / decade, logarithmic
Output impedance Load impedance Response time	2 x 4.7 Ohm, short circuit-proof > 10 kΩ < 30 ms
Transmitter identification	71.5 kΩ
HV transmitter alignment	at < 10 ⁻⁵ mbar
ATM transmitter alignment	at > 100 mbar
Semiconductor relay Adjustment range (N ₂) Hysteresis Switching characteristic Contact rating closed open Switching time	Switching thresholds SP 1/2 5.0 x 10 ⁻⁵ to 1500 mbar 10% of the threshold Low trip point < 30 V AC/DC, # 0.3 A resistive LED on LED off < 30 ms

WEBTank

Manuale d'USO

Codifica: MDU_WT-0519

Ed.01

Rev.04

06/03/21

Further Technical Data		THERMOVAC Transmitter TTR 101 / TTR 101 S2 / TTR 101 S2 PB	
Supply voltage at the transmitter	V DC	Class 2 / LPS +15 to +30	
Ripple		< 1 V _{pp}	
Power consumption	W	≤ 2,5	
Fuse to connect ahead		1 AT	
Electrical connection		FOC 68 socket	
Measurement cord		with shield: 0,14 mm ² / wire	
Cable length	m	≤ 100	
Earthing concept		--> "Electrical connection"	
Vacuum connection and signal earth		connected via 10 kV, 10 nF	
Materials exposed to vacuum			
Vacuum connection		Stainless steel 1.4435	
Heating filament		W	
Vacuum feedthrough		Glass	
Orifice (for DN 16 ISO-KF only)		Stainless steel	
Diaphragm		Ceramic	
Further materials		Ni, NiFe, Stainless steel 1.4301, SnAg	
Internal volume			
DN 16 ISO-KF	cm ³	4,7	
1/8- NPT	cm ³	5,2	
Admissible pressure	bar (abs.)	≤ 5	
Burst pressure	bar (abs.)	10	
Admissible temperature			
Operation (ambient)	°C	+10 to +50	
Bakeout at flange ^{*)} , max.	°C	≤ 80	
Heating filament	°C	< 180	
Storage	°C	-20 to +5	
Relative humidity			
Annual average		≤ 65% (non condensing)	
on 60 days		≤ 65% (non condensing)	
Mounting orientation		Horizontal and vertical	
Use		Indoors only, altitudes up to 2000 m (6562 ft.)	
Protection class	IP	40	
Weight	kg (lbs)	0,120 (0,265)	

*) For horizontal installation. While baking out, the technical data for measurement range, inaccuracy and repeatability may deviate

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
EU Declaration of Conformity

EU Declaration of Conformity

(Translation of original Declaration of Conformity)

The manufacturer: Leybold GmbH
 Bonner Strasse 498
 D-50968 Köln
 Germany

herewith declares that the products specified and listed below which we have placed on the market, comply with the applicable EU Directives. This declaration becomes invalid if modifications are made to the product without agreement of Leybold GmbH.

Product designation: THERMOVAC Transmitter

Type designation: TTR 101, TTR 101 S2, TTR 101 P8

Part numbers: 230350V01, 230351V01, 230352V01, 230353V01, 230354V01,
 230355V01, 230356V01, 230357V01, 230358V01

The products comply with the following Directives:



Electromagnetic Compatibility (2014/30/EU)

RoHS-Directive (2011/65/EU)

The following harmonized standards have been applied:

EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control, and laboratory use — Part 1: General requirements
EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use — EMC requirements — Part 1: General requirements Emissions: Group 1, Class B Immunity: Industrial electromagnetic environment
EN 60581:2012	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

Documentation office: Leybold GmbH, Bonner Straße 498, D-50968 Köln
 T: +49(0)221 347 0
 F: +49(0)221 347 1250
 documentation@leybold.com

<p>Cologne, August 30th, 2019</p>  <p>Th. Dreier Engineering Manager Product Company Cologne</p>	<p>Cologne, August 30th, 2019</p>  <p>Axel Guddes General Manager Product Company Cologne</p>
---	--

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Raspberry Pi

Raspberry Pi 2, Model B

Product Name	Raspberry Pi 2, Model B
Product Description	The Raspberry Pi 2 delivers 6 times the processing capacity of previous models. This second generation Raspberry Pi has an upgraded Broadcom BCM2836 processor, which is a powerful ARM Cortex-A7 based quad-core processor that runs at 900MHz. The board also features an increase in memory capacity to 1Gbyte.
RS Part Number	832-6274
Specifications	
Chip	Broadcom BCM2836 SoC
Core architecture	Quad-core ARM Cortex-A7
CPU	900 MHz
GPU	Dual Core VideoCore IV Multimedia Co-Processor Provides Open GL ES 2.0, hardware-accelerated OpenVG, and 1080p30 H.264 high-profile decode Capable of 1Gpixel/s, 1.5Gpixel/s or 24GROPs with texture filtering and DMA infrastructure
Memory	1GB LPDDR2
Operating System	Boots from Micro SD card, running a version of the Linux operating system.
Dimensions	85 x 56 x 17mm
Power	Micro USB socket 5V, 2A
Connectors:	
Ethernet	10/100 BaseT Ethernet socket
Video Output	HDMI (rev 1.3 & 1.4) Composite RCA (PAL and NTSC)
Audio Output	3.5mm jack, HDMI
USB	4 x USB 2.0 Connector
GPIO Connector	40-pin 2.54 mm (100 mil) expansion header: 2x20 strip Providing 27 GPIO pins as well as +3.3 V, +5 V and GND supply lines
Camera Connector	15-pin MIPI Camera Serial Interface (CSI-2)
JTAG	Not populated
Display Connector	Display Serial Interface (DSI) 15 way flat flex cable connector with two data lanes and a clock lane
Memory Card Slot	Micro SDIO

	<h1 style="text-align: center;">WEBTank</h1> <h2 style="text-align: center;">Manuale d'USO</h2> <p style="text-align: center;">Codifica: MDU_WT-0519</p>	Ed.01	Rev.04	06/03/21
---	--	-------	--------	----------

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Raspberry Pi Trading
30 Station Road
Cambridge
CB1 2JH

Web
<http://raspberrypi.org>

European Declaration of Conformity

Raspberry Pi (Trading) hereby declares under our own responsibility that the following product:

Raspberry Pi 2 Model B

Are in conformity with the following applicable community harmonised legislation:

Electromagnetic Compatibility Directive (EMC) 2014/30/EU,
Restriction of Hazardous Substances (RoHS) Directive 2011/65/EU

The following harmonised standards have been used to demonstrate conformity to these standards:

55032:2012 Class B
EN55024:2010,
EN61000-3-2:2014
EN61000-3-3:2013

Signed on behalf of the Raspberry Pi (Trading) Limited.


Gordon Hollingworth
Director of Engineering
Raspberry Pi (Trading) Limited

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High Efficiency 1.2MHz 2A Step Up Converter

FEATURES

- Integrated 80mΩ Power MOSFET
- 2V to 24V Input Voltage
- 1.2MHz Fixed Switching Frequency
- Internal 4A Switch Current Limit
- Adjustable Output Voltage
- Internal Compensation
- Up to 28V Output Voltage
- Automatic Pulse Frequency Modulation Mode at Light Loads
- up to 97% Efficiency
- Available in a 6-Pin SOT23-6 Package

APPLICATIONS

- Battery-Powered Equipment
- Set-Top Boxed
- LCD Bias Supply
- DSL and Cable Modems and Routers
- Networking cards powered from PCI or PCI express slots

GENERAL DESCRIPTION

The MT3608 is a constant frequency, 6-pin SOT23 current mode step-up converter intended for small, low power applications. The MT3608 switches at 1.2MHz and allows the use of tiny, low cost capacitors and inductors 2mm or less in height. Internal soft-start results in small inrush current and extends battery life.

The MT3608 features automatic shifting to pulse frequency modulation mode at light loads. The MT3608 includes under-voltage lockout, current limiting, and thermal overload protection to prevent damage in the event of an output overload. The MT3608 is available in a small 6-pin SOT-23 package.

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12W AC-DC Reliable Green Medical Adaptor

GSM12E series

(Cable connection)

(USB connection)

■ Features

- Universal AC input / Full range
- Medical safety approved (2 x MOPP) according to EN60601-1/EN60601-1-11
- Extremely low leakage current
- No load power consumption < 0.075W (< 0.1W for 18V/48V)
- **Energy efficiency Level VI**
- -20~+70°C wide range working temperature
- Class II power (no earth pin)
- Protections: Short circuit / Overload / Over voltage
- Fully enclosed plastic case
- 3 years warranty

■ Applications

- Blood glucose meter
- Blood pressure meter
- Nebulizer
- Inhaler
- Portable medical device
- Sleep apnea devices

■ Description

GSM12E is a highly reliable, 12W wall-mounted style single-output green medical adaptor series, which is compact and convenient for carry. This product is equipped with the standard 2-pin European plug. GSM12E is a class II power unit (no FG), accepting the input range from 80VAC to 264VAC that it can satisfy the demands for various types of medical electrical devices. The circuitry design meets the international medical standards (2*MOPP), having an ultra low leakage current (<100µA), fitting the medical devices in direct electrical contact with the patients.

With the working efficiency up to 87% and the extremely low no-load power consumption below 0.075W (< 0.1W for 18V/48V), GSM12E is compliant with EU ErP and meet CoC version 5 (Except for 18V/48V). The supreme feature allows the adaptor to save the energy when it is under either the operating mode or the standby mode. The entire series is approved for international safety regulations; moreover, it adopts the 94V-0 flame retardant plastic case that it can effectively prevent users from

■ Model Encoding

GSM12E 05 - P1J

DC plug type

Output voltage

European 2 pin AC plug

Rated wattage

Series name

P1J: Plug for standard model, 2.1ψ × 5.5ψ × 11mm, C+ tuning fork type

USB: USB(type A) for 5V model only

Other option available by customer requested (See page 5~6)

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<h1>Certificate</h1>		
Certificate no. TA 50371281 01		
License Holder: MEAN WELL Enterprises Co., Ltd. No. 28, Wuguan 3rd Rd., Wugu Dist., New Taipei City 248 Taiwan	Manufacturing Plant: ZTW1--50054701 001	
Test report no.: ZTW1- 50054701 001	Client Reference: ER-06022017-2Co	
Tested to: IEC 60601-1:2005+A1 EN 60601-1:2006+A11+A1+A12		
Certified Product: (SWITCHING POWER SUPPLY)	License Fee - Units	
Type Designation : GSM12Eyzzzzz (MW)	5	
y stands for : 05, 07, 09, 12, 15, 18, 24, 28 or 48	1	
z stands for : A to Z, 0 to 9, hyphen or blank	1	
Rated Voltage : AC 100-240V, 50/60Hz		
Rated Current : 0.4-0.2A		
Output Voltage : see Appendix		
Output Current : see Appendix		
max. Ambient Temperature: 40°C		
Protection Class : II		
This certificate is valid for maximum 5 years.		
7		
Licensed Test mark: IEC 60601-1 EN 60601-1 www.tuv.com ID: 141904408	TÜV Rheinland Taiwan Ltd Signature Dipl.-Ing. W. Hsu	Date of Issue (day/mo/yr) 13/02/2017
<small>TÜV Rheinland Taiwan Ltd., 11F., No. 75B, Sec. 4, Rode Road, Songshan Dist., Taipei City 105, Taiwan, R.O.C. - Tel: (886) 2172-7000 - Fax: (886) 2528-0018</small>		

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
Contenuto / Table Of Contents

0	Release Control Record.....	2
1	Equipment under test (EUT)	3
1.1	EUT Identification.....	3
1.2	Port identification.....	3
1.3	EUT classification.....	4
1.4	Additional information related to testing	4
1.5	Modifications incorporated in E.U.T.	4
1.6	Auxiliary equipment	4
1.7	Primary functions of the EUT	5
1.8	Performance of equipment under test	5
1.9	Performance criteria	6
2	Applications of test specifications 742 711 11	7
2.1	Immunity Test according to EN 61326-1.....	7
2.2	Emissions Test according to EN 61326-1.....	8
2.3	Deviations from Test Specification.....	8
3	Test specifications.....	9
4	Operating test modes and test conditions	9
5	Summary of test results	10
6	Photographic documentation	32

0 Release Control Record

Test Report Number	Reason of Change	Date of Issue
0141_21_EMCTR-0	Original release	06/03/2021

1 Equipment under test (EUT)

1.1 EUT Identification

Description	PRESSURE SENSOR
Model name or No.	WEBTANK
Part number / Serial No.	C21001
Manufacturer	Wireless Sensor Networks s.r.l.
Country of manufacturer	ITALY
Label	

1.2 Port identification

Port	Description	Type	Cable Length (m)	Type of connector
P1	Enclosure port	Plastic	---	screw
P2	AC power port	230V 50Hz	>3m	plug
P3	DC power port	Port not present	---	---
P4	I/O signal/ control ports (including functional earth lines)	Sensor signal	<2m	wired
P5	I/O signal/ control ports connected directly to power supply network	Port not present	---	---
P6	Wired Network ports	Port not present	---	---

1.3 EUT classification

EQUIPMENT CATEGORY	
<input checked="" type="checkbox"/>	Electrical measurement and test equipment
<input type="checkbox"/>	Electrical control equipment
<input type="checkbox"/>	Electrical laboratory equipment
INTENDED OPERATING EM ENVIRONMENT	
<input type="checkbox"/>	Residential, commercial and light-industrial locations
<input checked="" type="checkbox"/>	Industrial locations
<input type="checkbox"/>	Controlled electromagnetic laboratories or test and measurement areas
EMISSION CLASS (rif. EN 55011)	
<input checked="" type="checkbox"/>	A
<input type="checkbox"/>	B

1.4 Additional information related to testing

Intended Installation:	<input checked="" type="checkbox"/> Wall Mounted <input type="checkbox"/> Table Top <input type="checkbox"/> Handheld <input type="checkbox"/> Floor standing <input type="checkbox"/> Wall
Power Supply Requirement(s):	230V 50Hz

1.5 Modifications incorporated in E.U.T.

The following items are the modifications introduced in the equipment under test:

- **Added wurth ferrite 742 711 11 on internal power supply cable (2 loop of the wire around the ferrite core)**

1.6 Auxiliary equipment

- None

1.7 Primary functions of the EUT

The following table describes the primary functions and the representative parameter of the equipment under test according the manufacturer specifications:

Primary function	Representative parameter
Pressure Sensor	Measured Pressure

1.8 Performance of equipment under test

With reference to the above specified primary functions, the following table defines the acceptable level of the performance or permissible loss of function and the observation mode for each representative parameter of the equipment under test according to the technical instructions by the manufacturer.

Representative parameter	Acceptable level of performance	Observation mode		
		Acquisition	Test equipment	Test n.
Measured Pressure	Tolerance of $\pm 50\%$ as declared by the manufacturer	Operator	Visual inspection	All immunity test

1.9 Performance criteria

6.4.1 General

The general principles (performance criteria) for the evaluation of the immunity test results are the following.

6.4.2 Performance criterion A

The equipment shall continue to operate as intended during and after the test. No degradation of performance or loss of function is allowed below a performance level specified by the manufacturer, when the equipment is used as intended. The performance level may be replaced by a permissible loss of performance. If the minimum performance level or the permissible performance loss is not specified by the manufacturer, either of these may be derived from the product description and documentation and what the user may reasonably expect from the equipment if used as intended.

6.4.3 Performance criterion B

The equipment shall continue to operate as intended after the test. No degradation of performance or loss of function is allowed below a performance level specified by the manufacturer, when the equipment is used as intended. The performance level may be replaced by a permissible loss of performance. During the test, degradation of performance is however allowed. No change of actual operating state or stored data is allowed. If the minimum performance level or the permissible performance loss is not specified by the manufacturer, either of these may be derived from the product description and documentation and what the user may reasonably expect from the equipment if used as intended.

EXAMPLE 1 A data transfer is controlled/checked by parity check or by other means. In the case of malfunctioning, such as caused by a lightning strike, the data transfer will be repeated automatically. The reduced data transfer rate at this time is acceptable.

EXAMPLE 2 During testing, an analogue function value may deviate. After the test, the deviation vanishes.

EXAMPLE 3 In the case of a monitor used only for man-machine monitoring, it is acceptable that some degradation takes place for a short time, such as flashes during the burst application.

EXAMPLE 4 An intended change of the operating state is allowed if self-recoverable.

6.4.4 Performance criterion C

Temporary loss of function is allowed, provided the function is self-recoverable or can be restored by the operation of the controls.

EXAMPLE 1 In the case of an interruption in the mains longer than the specified buffer time, the power supply unit of the equipment is switched off. The switch-on may be automatic or carried out by the operator.

EXAMPLE 2 After a programme interruption caused by a disturbance, the processor functions of the equipment stops at a defined position and is not left in a "crashed state". The operator's decision prompts may be necessary.

EXAMPLE 3 The test results in an opening of an over-current protection device that is replaced or reset by the operator.

2 Applications of test specifications

2.1 Immunity Test according to EN 61326-1

Table 1 – Immunity test requirements for equipment intended to be used in a basic electromagnetic environment				
Port	Phenomenon	Basic standard	Test value	Performance criterion
Enclosure	Electrostatic discharge (ESD)	IEC 61000-4-2	4 kV contact discharge 8 kV air discharge	B B
	Electromagnetic field	IEC 61000-4-3	3 V/m (80 MHz to 1 GHz) 3 V/m (1,4 GHz to 2 GHz) 1 V/m (2,0 GHz to 2,7 GHz)	A A A
	Power frequency magnetic field	IEC 61000-4-8	3 A/m (50 Hz, 80 Hz) ¹	A
AC power (including protective earth)	Voltage dip	IEC 61000-4-11	0 % during half cycle 0 % during 1 cycle 70 % during 25/30 ^c cycles	B B C
	Short interruptions	IEC 61000-4-11	0 % during 250/300 ^d cycles	C
	Burst	IEC 61000-4-4	1 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	0,5 kV ^b /1 kV ^b	B
	Conducted RF	IEC 61000-4-6	3 V (150 kHz to 80 MHz)	A
DC power ^e (including protective earth)	Burst	IEC 61000-4-4	1 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	0,5 kV ^b /1 kV ^b	B
	Conducted RF	IEC 61000-4-6	3 V (150 kHz to 80 MHz)	A
I/O signal/control (including functional earth)	Burst	IEC 61000-4-4	0,5 kV ^b (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	1 kV ^{b, c}	B
	Conducted RF	IEC 61000-4-6	3 V ^f (150 kHz to 80 MHz)	A
I/O signal/control connected directly to mains supply	Burst	IEC 61000-4-4	1 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	0,5 kV ^b /1 kV ^b	B
	Conducted RF	IEC 61000-4-6	3 V (150 kHz to 80 MHz)	A

^a Line to line.
^b Line to ground.
^c Only in the case of long-distance lines (see 3.10).
^d Only in the case of lines > 3 m.
^e For example "25/30 cycles" means "25 cycles for 50 Hz test" or "30 cycles for 60 Hz test".
^f Only to magnetically sensitive equipment. CRT display interference is allowed above 1 A/m.
^g DC connections between parts of equipment/system which are not connected to a d.c. distribution network are treated as I/O signal/control ports.

Table 2 – Immunity test requirements for equipment intended to be used in an industrial electromagnetic environment				
Port	Phenomenon	Basic standard	Test value	Performance criterion
Enclosure	Electrostatic discharge (ESD)	IEC 61000-4-2	4 kV contact discharge 8 kV air discharge	B B
	Electromagnetic field	IEC 61000-4-3	10 V/m (80 MHz to 1 GHz) 3 V/m (1,4 GHz to 2 GHz) 1 V/m (2,0 GHz to 2,7 GHz)	A A A
	Power frequency magnetic field	IEC 61000-4-8	30 A/m (50 Hz, 80 Hz) [*]	A
AC power (including protective earth)	Voltage dip	IEC 61000-4-11	0 % during 1 cycle 40 % during 10/12 ^g cycles 70 % during 25/30 ^d cycles	B C C
	Short interruptions	IEC 61000-4-11	0 % during 250/300 ^d cycles	C
	Burst	IEC 61000-4-4	2 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	1 kV ^b /2 kV ^b	B
	Conducted RF	IEC 61000-4-6	3 V ^f (150 kHz to 80 MHz)	A
DC power ^g (including protective earth)	Burst	IEC 61000-4-4	2 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	1 kV ^b /2 kV ^b	B
	Conducted RF	IEC 61000-4-6	3 V ^f (150 kHz to 80 MHz)	A
I/O signal/control (including functional earth)	Burst	IEC 61000-4-4	1 kV (5/50 ns, 5 kHz) ^f	B
	Surge	IEC 61000-4-5	1 kV ^{b, c}	B
	Conducted RF	IEC 61000-4-6	3 V ^f (150 kHz to 80 MHz)	A
I/O signal/control connected directly to mains supply	Burst	IEC 61000-4-4	2 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	1 kV ^b /2 kV ^b	B
	Conducted RF	IEC 61000-4-6	3 V ^f (150 kHz to 80 MHz)	A

^a Line to line.
^b Line to ground.
^c Only in the case of long-distance lines (see 3.10).
^d Only in the case of lines > 3 m.
^e Only to magnetically sensitive equipment. CRT display interference is allowed above 1 A/m.
^f DC connections between parts of equipment/system which are not connected to a d.c. distribution network are treated as I/O signal/control ports.
^g For example "25/30 cycles" means "25 cycles for 50 Hz test" or "30 cycles for 60 Hz test".

Table 3 – Immunity test requirements for equipment intended to be used in a controlled electromagnetic environment				
Port	Phenomenon	Basic standard	Test value	Performance criterion
Enclosure	Electrostatic discharge (ESD)	IEC 61000-4-2	4 kV contact discharge 8 kV air discharge	B B
	Electromagnetic field	IEC 61000-4-3	1 V/m (80 MHz to 1 GHz) 1 V/m (1,4 GHz to 2 GHz) 1 V/m (2,0 GHz to 2,7 GHz)	A A A
AC power (including protective earth)	Voltage dip	IEC 61000-4-11	0 % during half cycle	B
	Burst	IEC 61000-4-4	1 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	0,5 kV ^b /1 kV ^b	B
	Conducted RF	IEC 61000-4-6	1 V (150 kHz to 80 MHz)	A
DC power ^a (including protective earth)	Burst	IEC 61000-4-4	1 kV (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	Not required	-
	Conducted RF	IEC 61000-4-6	1 V (150 kHz to 80 MHz)	A
I/O signal/control (including functional earth)	Burst	IEC 61000-4-4	0,5 kV ^b (5/50 ns, 5 kHz)	B
	Surge	IEC 61000-4-5	Not required	-
	Conducted RF	IEC 61000-4-6	1 V ^f (150 kHz to 80 MHz)	A

^a Line to line.
^b Line to ground.
^c Only in the case of lines > 3 m.
^d DC connections between parts of equipment/system which are not connected to a d.c. distribution network are treated as I/O signal/control ports.

--- BLANK ---

2.2 Emissions Test according to EN 61326-1

7.2 Emission limits

The equipment shall be classified and respective information shall be provided per the applicable group and class as specified within CISPR 11:2009, Clause 5. Equipment classification and choice of respective limits shall be determined after taking into account the intended environment and emission requirement in the areas of use.

For Class A equipment, the limits, the measuring methods and the provisions given in CISPR 11 apply.

For Class B equipment, the limits, the measuring methods and the provisions given in CISPR 11, [EC 61000-3-2](#) (or [EC 61000-3-12](#)) and [EC 61000-3-3](#) (or [EC 61000-3-11](#)) apply.

For equipment using frequencies in the ISM bands, see CISPR 11.

2.3 Deviations from Test Specification

For the measurements contained within this test report, there were no deviations from, additions to, or exclusions from the test specification identified above, nor from the requirements defined in the basic standards called up within it.

3 Test specifications

REFERENCE EUROPEAN STANDARDS	
EN 61326-1:2013	Electrical equipment for measurement, control and laboratory use – EMC requirements - Part 1: General requirements

4 Operating test modes and test conditions

In the following table there are the operating conditions adopted during tests identified by an indicator (#..) at which has been referred the item “Operating condition of the equipment under test”.

Operating condition	Description
#1	CONTINUOUS PRESSURE MEASUREMENT

5 Summary of test results

5.1 Emission tests

Port		Test	Reference standard	Operating condition	Results
P1	Enclosure	Radiated emission	CISPR 11	#1	Within the limits
P2	AC power port	Harmonic current emissions	IEC 61000-3-2	#1	Within the limits
		Voltage fluctuations flicker	IEC 61000-3-3	#1	Within the limits
		Conducted emission	CISPR 11 CISPR 22	#1	Within the limits

5.2 Immunity tests

Port		Phenomena	Reference standard	Operating condition	Result
1	Enclosure	RF electromagnetic field	EN 61000-4-3	#1	COMPLIANT
		Electrostatic Discharge (ESD)	EN 61000-4-2	#1	COMPLIANT
		Power frequency magnetic field	EN 61000-4-8	---	Not Applicable ¹
2	AC MAINS ports	Fast transients	EN 61000-4-4	#1	COMPLIANT
		RF common mode	EN 61000-4-6	#1	COMPLIANT
		Surge	EN 61000-4-5	#1	COMPLIANT
		Voltage Dips / Interruptions	EN 61000-4-11	#1	COMPLIANT
3	Input and output DC POWER ports	Fast transients	EN 61000-4-4	---	Not Applicable ²
		RF common mode	EN 61000-4-6		
		Surge	EN 61000-4-5		
4	I/O signal/ control ports (including functional earth lines)	Fast transients	EN 61000-4-4	---	Not Applicable ²
		RF common mode	EN 61000-4-6		
		Surge	EN 61000-4-5		
5	I/O signal/ control ports connected directly to power supply network	Fast transients	EN 61000-4-4	---	Not Applicable ²
		RF common mode	EN 61000-4-6		
		Surge	EN 61000-4-5		
6	Wired Network Ports	Fast transients	EN 61000-4-4	---	Not Applicable ²
		RF common mode	EN 61000-4-6		

¹ The EUT doesn't contain significant elements that are susceptible to magnetic fields

² Port not present

TEST 1

Voltage Disturbances

Reference Document: CISPR 11; CISPR 22

Test Setup:

Test Location:

Test Equipment used for test:

Measurement Uncertainty

In acc. to ref. doc.

RadioMotive s.r.l.

EMC.023 LISN R&S ESH2-Z5

EMC.018 EMI RECEIVER R&S ESW8

2,681dB


Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1


Test results: **WITHIN THE LIMITS**

NOTE: EUT is a class A equipment but conducted emissions respects also class B limits

LINE: PHASE (L1)

LINE: NEUTRAL (N)

TEST 2

Radiated electromagnetic disturbances

Reference Document: CISPR 11

Test Setup:

In acc. to ref. doc.

Test Location:

RadioMotive s.r.l.

Test Equipment used for test:

EMC.048 Antenna AAronia mod. Bicolog 30100E

EMC.037 Miteq Preamplifier

EMC.018 EMI RECEIVER R&S ESW8


Measurement Uncertainty

4,82dB

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %Rh	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1

Test results: **WITHIN THE LIMITS**

Final Result

Frequency (MHz)	QuasiPeak (dBµV/m)	Limit (dBµV/m)	Margin (dB)	Height (cm)	Pol	Azimuth (deg)
118.530000	46.81	50.00	3.19	155.0	H	180.0
177.510000	45.82	50.00	4.18	155.0	H	180.0
187.680000	45.53	50.00	4.47	155.0	H	180.0

TEST 3

Harmonics Current Emissions

Reference Document: IEC 61000-3-2

Test Setup:

Test Location:

Test Equipment used for test:

In acc. to ref. doc.

RadioMotive s.r.l.

EMC.060 Digital Power Analyzer DPA 500N

EMC.059 Power Supply Zenone GTS300GL

Measurement Uncertainty

2,681dB

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V ~ 50Hz

Operating Condition (Rif. Section 4): #1

Test results: **WITHIN THE LIMITS**

Measure Results	
Standard Specific Results for IEC 61000-3-2 (Edition 5)	
Standard Group:	Industry
Standard Name:	IEC 61000-3-2 (Edition 5) Limits for harmonic current emissions (equipment input current < 16 A per phase)
Device Under Test:	PASS
Power Source:	PASS
Connection Type:	L - N
Classification:	Class A
Appl. of Limits:	less than or equal to 150 % (Without POHC Enhancement)
Current limits are disabled because rated power is less than 75W.	
Check Harmonics 2..40	
<i>First detected harmonic order > 150 %</i>	
Line 1:	None
<i>Harmonics orders > 150 %</i>	
Line 1:	None
<i>Harmonics orders with average > 100 %</i>	
Line 1:	None
Measured values:	
<i>Fundamental Current</i>	
Line 1:	0,024 A
<i>Active input Power</i>	
Line 1:	5,4 W *
<i>Circuit power factor</i>	
Line 1:	0,387 *
* Absolute value.	

For equipment with a total rated power < 75W limits are not defined in the standard.

TEST 4

Limit Voltage Fluctuations and Flicker Emissions

Reference Document: IEC 61000-3-3

Test Setup:

Test Location:

Test Equipment used for test:

Measurement Uncertainty

In acc. to ref. doc.

RadioMotive s.r.l.

Digital Power Analyzer DPA 500N

Elettrotec TPS/M

2,681dB

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V ~ 50Hz

Operating Condition (Rif. Section 4): #1

Test results: WITHIN THE LIMITS

Flicker Results									
Standard Specific Results for IEC 61000-3-3 (Edition 3)									
Standard Group:	Industry								
Standard Name:	IEC 61000-3-3 (Edition 3) Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection								
Test Condition:	General Test Conditions								
Analysis Status:	PASS								
Flicker Measurements Settings									
Main line:	230V, 50Hz								
Flicker Meter:	230V / 50Hz								
Flicker Impedance:	Zref								
Observation Time:	1 x 10 min								
Measurements done:	1								
Flicker Measurements									
	Pst	Max Pst	Max Dc	Max Dmax	Max Tmax				
Line 1:	0,012	0,028	0	0,246	0				
Limits:	0,65	1	3,3	4	0,5				
Results:	PASS	PASS	PASS	PASS	PASS				
Pst Data									
Flicker (Line 1)									
Meas. No.	P0,1	P1s	P3s	P10s	P50s	Pst	dc [%]	dmax [%]	Tmax [s]
1	0	0	0	0,001	0,005	0,028	0	0,246	0
Short-term Flicker Severity (Pst) (Line 1)									
									
Long-term Flicker Severity (Pft) (Line 1)									
									

TEST 5

Immunity to Electrostatic Discharge

Reference Document: IEC 61000-4-2

Test Setup:

In acc. to ref. doc.

Test Location:

RadioMotive s.r.l.

Test Equipment used for test:

EMC.009 EMtest ESD Generator - DITO

Measurement Uncertainty

VOLTAGE DISCHARGE 3.69%

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1

Test results: Compliant

CONTACT DISCHARGE (for conductive surfaces)

- Level: ± 4 kV
- Number of discharge: 10 positive, 10 negative
- Repetition rate: 1Hz

TEST RESULTS		PERFORMANCE CRITERIA			RESULT
Discharge point	Pol.	COMPLIANT		NOT COMPLIANT	
		A	B	C	
Metallic plate on upper side	+	10	---	---	No performance degradation
	-	10			

N: No discharge occurs

AIR DISCHARGE (for not conductive surfaces)

- Level : ± 8 kV
- Number of discharge: 10 positive, 10 negative
- Repetition rate: 1Hz

TEST RESULTS		PERFORMANCE CRITERIA			RESULT
Discharge point	Pol.	COMPLIANT		NOT COMPLIANT	
		A	B	C	
Enclosure all sides	+	N			No performance degradation
	-	N			

N: No discharge occurs

Indirect DISCHARGE to VCP

- Level : ± 4 kV
- Number of discharge: 10 positive, 10 negative
- Repetition rate: 1Hz

TEST RESULTS		PERFORMANCE CRITERIA			RESULT
Discharge point	Pol.	COMPLIANT		NOT COMPLIANT	
		A	B	C	
VCP to Enclosure sides	+	N			No performance degradation
	-	N			

Indirect DISCHARGE to hCP

- Level : ± 4 kV
- Number of discharge: 10 positive, 10 negative
- Repetition rate: 1Hz

TEST RESULTS		PERFORMANCE CRITERIA			RESULT
Discharge point	Pol.	COMPLIANT		NOT COMPLIANT	
		A	B	C	
VCP to Enclosure bottom side	+	N			No performance degradation
	-	N			

TEST 6

Immunity to Radiated RF Electromagnetic Fields

Reference Document: IEC 61000-4-3

Test Setup:	In acc. to ref. doc.
Test Location:	RadioMotive s.r.l.
Test Equipment used for test:	EMC.013 Signal Generator – HP 8673D EMC.033 RF Amplifier KALMUS 7200LC-/1-70-826-001 EMC.035 RF Amplifier BONN BLMA 1060-30 EMC.025 Power Meter Agilent E4419B EMC.026 Power Sensor E4412A EMC.027 Power Sensor E4412A EMC.030 Antenna Log Periodic PROTEL ARL70F2500XZ EMC.028 Antenna Horn RF-SPIN DRH18-EX
Measurement Uncertainty	Level 23.15%

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50 Hz

Operating Condition (Rif. Section 4): #1

Test results: **COMPLIANT**

TEST RESULTS

80-2700 MHz AM, 80% 1kHz, 10V/m

POLAR. PORT	VERTICAL			HORIZONTAL			RESULT
	COMPLIANT A	NOT COMPLIANT B C		COMPLIANT A	NOT COMPLIANT B C		
ENCLOSURE front side	√			√			No performance degradation
ENCLOSURE rear side	√			√			
ENCLOSURE right side	√			√			
ENCLOSURE left side	√			√			

TEST 7

Immunity to Fast Transients (Bursts)

Reference Document: IEC 61000-4-4


Test Setup:	In acc. to ref. doc.
Test Location:	RadioMotive s.r.l.
Test Equipment used for test:	EMC.010 Burst generator EMTest compact NX5 bsp-1-300-16
Measurement Uncertainty	VOLTAGE 3.82%

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1

Test results: **COMPLIANT**

AC MAINS

Test Parameters	Values
Voltage (V)	1000 V
Polarity (Pol)	Alternate
Repetition frequency (f)	5 kHz
Burst duration (td)	15 ms
Burst period (tr)	300 ms
Duration (T)	60 s
Coupling (Cpl)	L, N, PE, L+N, L+PE, N+PE, L+N+PE
Phase synchronization (Sync)	Off
Graphic representation of burst	 <p>The diagram shows a voltage waveform V over time. It consists of two bursts of high-frequency pulses. Each burst is followed by a period of zero voltage. The duration of each burst is labeled as 'td' and the duration of the zero-voltage period is labeled as 'tr'. The total duration of one burst and its following zero-voltage period is labeled as 'T'.</p>
Test Results	
Result	PASS
Note	No performance degradation

TEST 8

Immunity to Surge

Reference Document: IEC 61000-4-5

Test Setup:	In acc. to ref. doc.
Test Location:	RadioMotive s.r.l.
Test Equipment used for test:	EMC.010 Surge generator EMTest compact NX5 bsp-1-300-16
Measurement Uncertainty	Voltage 7.33%

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1

Test results: **COMPLIANT**

TEST RESULTS

AC MAINS

TEST	COUPLING MODE	TEST VOLTAGE	TEST FREQUENCY	NUMBER OF PULSES	PHASE ANGLE (φ)
1	L1 - PE (COMMON MODE)	+/-2kV	1pul./min.	5	0°, 90°, 180°, 270°
2	N - PE (COMMON MODE)	+/-2kV	1pul./min.	5	0°, 90°, 180°, 270°
3	L1 - N (DIFF. MODE)	+/-1kV	1pul./min.	5	0°, 90°, 180°, 270°
RESULT	Performance Criteria			NOTES	
	COMPLIANT		NOT COMPLIANT		
	A	B	C		
TEST 1	✓			No performance degradation	
TEST 2	✓			No performance degradation	
TEST 3	✓			No performance degradation	

TEST 9

Immunity to Radio-Frequency common mode disturbances

Reference Document: IEC 61000-4-6

Test Setup:

In acc. to ref. doc.

Test Location:

RadioMotive s.r.l.

Test Equipment used for test:

EMC.011 Test system Teseq NSG 4070B-80

EMC.006 CDN Teseq M016

Measurement Uncertainty

1.55%

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1

Test results: **COMPLIANT**

TEST RESULTS

AM 80% 1kHz 3V

PORT n.		PERFORMANCE CRITERIA			RESULT
		COMPLIANT	NOT COMPLIANT		
		A	B	C	
P2	AC power ports	✓			No performance degradation
P3	DC power ports	---			Not applicable: port not present
P4	I/O signal/ control (including functional earth)	---			Not applicable: Cable length <3m
P5	I/O signal/ control connected directly to mains supply	---			Not applicable: port not present

TEST 10

Immunity to Voltage Dips / Short Interruptions

Reference Document: EN 61000-4-11

Test Setup:	In acc. to ref. doc.
Test Location:	RadioMotive s.r.l.
Test Equipment used for test:	EMC.010 Power fail generator EMTest compact NX5 bsp-1-300-16
Measurement Uncertainty	Voltage 3.86%

Test conditions:	Measured
Ambient temperature: 23°C ±5°C	24 °C
Ambient humidity: 25 – 75 %rH	45 %
Pressure: 85 – 106 kPa – (860 mbar – 1060 mbar)	960 mbar
Voltage	230V 50Hz

Operating Condition (Rif. Section 4): #1

Test results: **COMPLIANT**

TEST RESULTS

Voltage Dips

TEST LEVEL			COMPLIANT		NOT COMPLIANT	RESULT
% res.	cycles	φ (°)	A	B	C	
0 %	1	0° 90° 180° 270°	✓ ✓ ✓ ✓			No performance degradation

TEST LEVEL			COMPLIANT		NOT COMPLIANT	RESULT
% res.	cycles	φ (°)	A	B	C	
40 %	10	0° 90° 180° 270°	✓ ✓ ✓ ✓			No performance degradation

TEST LEVEL			COMPLIANT		NOT COMPLIANT	RESULT
% res.	cycles	φ (°)	A	B	C	
70 %	25	0° 90° 180° 270°	✓ ✓ ✓ ✓			No performance degradation

Voltage interruption

TEST LEVEL			COMPLIANT			NOT COMPLIANT	RESULT
% res.	cycles	φ (°)	A	B	C		
0 %	250	0° 90° 180° 270°			✓ ✓ ✓ ✓	EUT shutdown. Automatic restart of operating cycle	

6 Photographic documentation

PHOTO 1 – EUT IDENTIFICATION

PHOTO 2 – RADIATED EMISSIONS DISTURBANCES

PHOTO 3 – CONDUCTED EMISSIONS DISTURBANCES

PHOTO 4 – RADIATED IMMUNITY SETUP


PHOTO 5 – BURST/SURGE/DIPS IMMUNITY SETUP

PHOTO 6 – RF CM IMMUNITY SETUP


PHOTO 7 – BURST/SURGE/DIPS IMMUNITY SETUP

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Via Tevere 63 - 22073

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Tel. +39 031 4134033

info@radiomotive.it
www.radiomotive.it

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12 300544666-001, Rev D - 2019 - © Leybold
RadioMotive Srl